

Cleaved-Caspase-7 (S199) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10005	SYNONYMS
SIZE	100 ug	CASP7; MCH3; Caspase-7; CASP-7; Apoptotic protease
HOST	Rabbit	Mch-3; CMH-1; ICE-like apoptotic protease 3; ICE-LAP3

PRODUCT INFORMATION

GENE NAME	CASP7
PROTEIN NAME	Caspase-7
IMMUNOGEN	Synthesized peptide derived from the Internal region of human Caspase-7. at AA range: 150-230

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.