

Cleaved-Caspase-9 (D353) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10007	SYNONYMS
SIZE	100 ug	CASP9; MCH6; Caspase-9; CASP-9; Apoptotic protease
HOST	Rabbit	Mch-6; Apoptotic protease-activating factor 3; APAF-3; ICE-like apoptotic protease 6; ICE-LAP6

PRODUCT INFORMATION

GENE NAME	CASP9
PROTEIN NAME	Caspase-9
IMMUNOGEN	Synthesized peptide derived from the C-terminal region of human Caspase-9. at AA range: 280-360

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Mouse,Rat
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.