

His-Tag Mouse Monoclonal Antibody

OVERVIEW

CATALOG NO.	PTX1001
SIZE	100 µL
HOST	Mouse
ISOTYPE	IgG1
CLONALITY	Monoclonal Antibody

APPLICATIONS AND STORAGE

APPLICATIONS	WB, IP, IF
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles
DELIVERY CONDITION	Blue ice