

Adiponectin Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10010	SYNONYMS
SIZE	100 ug	ADIPOQ; ACDC; ACRP30; APM1; GBP28;
HOST	Rabbit	Adiponectin; 30 kDa adipocyte complement-related protein; Adipocyte complement-related 30 kDa protein; ACRP30; Adipocyte; C1q and collagen domain-containing protein; Adipose most abundant gene transcript 1

PRODUCT INFORMATION

GENE NAME	ADIPOQ
PROTEIN NAME	Adiponectin
IMMUNOGEN	Synthesized peptide derived from the N-terminal region of human Adiponectin. at AA range: 10-90

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.