

Cleaved-Caspase-2 p18 (G170) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10015	SYNONYMS
SIZE	100 ug	CASP2; ICH1; NEDD2; Caspase-2; CASP-2; Neural precursor cell expressed developmentally down-regulated protein 2; NEDD-2; Protease ICH-1
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	CASP2
PROTEIN NAME	Caspase-2
IMMUNOGEN	Synthesized peptide derived from the Internal region of human Caspase-2 p18. at AA range: 120-200

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.