

Cleaved-Cathepsin C HC (R394) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10024	SYNONYMS
SIZE	100 ug	CTSC; CPPI; Dipeptidyl peptidase 1; Cathepsin C;
HOST	Rabbit	Cathepsin J; Dipeptidyl peptidase I; DPP-I; DPPI; Dipeptidyl transferase

PRODUCT INFORMATION

GENE NAME	CTSC
PROTEIN NAME	Dipeptidyl peptidase 1
IMMUNOGEN	Synthesized peptide derived from the C-terminal region of human Cathepsin C HC. at AA range: 320-400

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.