

Cleaved-Cathepsin D LC (G65) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10026	SYNONYMS
SIZE	100 ug	CTSD; CPSD; Cathepsin D
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	CTSD
PROTEIN NAME	Cathepsin D
IMMUNOGEN	Synthesized peptide derived from the Internal region of human Cathepsin D LC. at AA range: 20-100

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Monkey
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.