

Cleaved-Cathepsin L2 (L114) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10031	SYNONYMS
SIZE	100 ug	CTSL2; CATL2; CTSU; CTSV; Cathepsin L2;
HOST	Rabbit	Cathepsin U; Cathepsin V

PRODUCT INFORMATION

GENE NAME	CTSL2
PROTEIN NAME	Cathepsin L2
IMMUNOGEN	Synthesized peptide derived from the Internal region of human Cathepsin L2. at AA range: 70-150

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Monkey
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.