

Cleaved-Cathepsin Z (L62) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10032	SYNONYMS
SIZE	100 ug	CTSZ; Cathepsin Z; Cathepsin P; Cathepsin X
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	CTSZ
PROTEIN NAME	Cathepsin Z
IMMUNOGEN	Synthesized peptide derived from the Internal region of human Cathepsin Z. at AA range: 20-100

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Monkey
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.