

## Cleaved-Lamin A (D230) Polyclonal Antibody

### OVERVIEW

---

<b>CATALOG NO.</b>	PTX10040	<b>SYNONYMS</b>
<b>SIZE</b>	100 ug	LMNA; LMN1; Prelamin-A/C
<b>HOST</b>	Rabbit	

### PRODUCT INFORMATION

---

<b>GENE NAME</b>	LMNA
<b>PROTEIN NAME</b>	Prelamin-A/C
<b>IMMUNOGEN</b>	Synthesized peptide derived from the Internal region of human Lamin A. at AA range: 150-230

### APPLICATIONS AND STORAGE

---

<b>SPECIES REACTIVITY</b>	Human,Mouse,Rat
<b>TESTED APPLICATIONS</b>	WB,IHC,ELISA
<b>RECOMMENDED DILUTION</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
<b>DELIVERY CONDITIONS</b>	Blue Ice
<b>STABILITY &amp; STORAGE</b>	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

**For research use only.**