

Cleaved-MMP-14 (Y112) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10043	SYNONYMS
SIZE	100 ug	MMP14; Matrix metalloproteinase-14; MMP-14; MMP-X1; Membrane-type matrix metalloproteinase 1; MT-MMP 1; MTMMP1; Membrane-type-1 matrix metalloproteinase; MT1-MMP; MT1MMP
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	MMP14
PROTEIN NAME	Matrix metalloproteinase-14
IMMUNOGEN	Synthesized peptide derived from the Internal region of human MMP-14. at AA range: 70-150

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.