

Cleaved-MMP-15 (Y132) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10044	SYNONYMS
SIZE	100 ug	MMP15; Matrix metalloproteinase-15; MMP-15;
HOST	Rabbit	Membrane-type matrix metalloproteinase 2; MT-MMP 2; MTMMP2; Membrane-type-2 matrix metalloproteinase; MT2-MMP; MT2MMP; SMCP-2

PRODUCT INFORMATION

GENE NAME	MMP15
PROTEIN NAME	Matrix metalloproteinase-15
IMMUNOGEN	Synthesized peptide derived from the Internal region of human MMP-15. at AA range: 90-170

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Monkey
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.