

Cleaved-Notch 1 (V1754) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10049	SYNONYMS
SIZE	100 ug	NOTCH1; TAN1; Neurogenic locus notch homolog
HOST	Rabbit	protein 1; Notch 1; hN1; Translocation-associated notch protein TAN-1

PRODUCT INFORMATION

GENE NAME	NOTCH1
PROTEIN NAME	Neurogenic locus notch homolog protein 1
IMMUNOGEN	Synthesized peptide derived from the Internal region of human Notch 1. at AA range: 1710-1790

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.