

Cleaved-Notch 2 (D1733) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10050	SYNONYMS
SIZE	100 ug	NOTCH2; Neurogenic locus notch homolog protein 2;
HOST	Rabbit	Notch 2; hN2

PRODUCT INFORMATION

GENE NAME	NOTCH2
PROTEIN NAME	Neurogenic locus notch homolog protein 2
IMMUNOGEN	Synthesized peptide derived from the Internal region of human Notch 2. at AA range: 1660-1740

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.