

Cleaved-SUMO-2/3 (G93) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10053	SYNONYMS
SIZE	100 ug	SUMO2; SMT3A; SMT3H2; Small ubiquitin-related
HOST	Rabbit	modifier 2; SUMO-2; HSMT3; SMT3 homolog 2; SUMO-3; Sentrin-2; Ubiquitin-like protein SMT3A; Smt3A

PRODUCT INFORMATION

GENE NAME	SUMO2
PROTEIN NAME	Small ubiquitin-related modifier 2
IMMUNOGEN	Synthesized peptide derived from the Internal region of human SUMO-2/3. at AA range: 20-100

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.