

Cleaved-Tumstatin (P1426) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10055	SYNONYMS
SIZE	100 ug	COL4A3; Collagen alpha-3(IV) chain; Goodpasture antigen
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	COL4A3
PROTEIN NAME	Collagen alpha-3(IV) chain
IMMUNOGEN	Synthesized peptide derived from the C-terminal region of human Tumstatin. at AA range: 1380-1460

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Monkey
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.