

HABP2 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10060	SYNONYMS
SIZE	100 ug	HABP2; HGFAL; PHBP; Hyaluronan-binding protein 2;
HOST	Rabbit	Factor VII-activating protease; Factor seven-activating protease; FSAP; Hepatocyte growth factor activator-like protein; Plasma hyaluronan-binding protein

PRODUCT INFORMATION

GENE NAME	HABP2
PROTEIN NAME	Hyaluronan-binding protein 2
IMMUNOGEN	Synthesized peptide derived from the Internal region of human HABP2. at AA range: 270-350

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.