

Cleaved-Integrin $\alpha 6$ LC (E942) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10062	SYNONYMS
SIZE	100 ug	ITGA6; Integrin alpha-6; CD49 antigen-like family member F; VLA-6; CD antigen CD49f
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	ITGA6
PROTEIN NAME	Integrin alpha-6
IMMUNOGEN	Synthesized peptide derived from the C-terminal region of human Integrin $\alpha 6$ LC. at AA range: 900-980

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.