

Cleaved-Plasminogen HC A short form (V98) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10072	SYNONYMS
SIZE	100 ug	PLG; Plasminogen
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	PLG
PROTEIN NAME	Plasminogen
IMMUNOGEN	Synthesized peptide derived from the N-terminal region of human Plasminogen HC A short form. at AA range: 50-130

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.