

Mono-Methyl-Histone H3 (K37) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10079	SYNONYMS
SIZE	100 ug	HIST1H3A; H3FA; HIST1H3B; H3FL; HIST1H3C;
HOST	Rabbit	H3FC; HIST1H3D; H3FB; HIST1H3E; H3FD; HIST1H3F; H3FI; HIST1H3G; H3FH; HIST1H3H; H3FK; HIST1H3I; H3FF; HIST1H3J; H3FJ; Histone H3.1; Histone H3.1; Histone H3.2; Histone H3/m; Histone H3/o; H3F3A; H3.3A; H3F3; PP781;

PRODUCT INFORMATION

GENE	HIST1H3A/HIST1H3/HIST1H3C/HIST1H3D/HIST1H3E/HIST1H3F/HIST1H3G/HIST1H3H/HIST1H3I/H
NAME	IST1H3J/HIST2H3A/HIST2H3C/HIST2H3D/H3F3A/H3F3B/H3F3C
PROTEI	Histone H3.1/Histone H3.2/Histone H3.3/Histone H3.3C
N NAME	
IMMUN	Synthesized peptide derived from human Histone H3 around the mono-methylation site of K37.
OGEN	

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.