

Acetyl-NFκB-p65 (K310) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10093	SYNONYMS
SIZE	100 ug	RELA; NFKB3; Transcription factor p65; Nuclear factor
HOST	Rabbit	NF-kappa-B p65 subunit; Nuclear factor of kappa light polypeptide gene enhancer in B-cells 3

PRODUCT INFORMATION

GENE NAME	RELA
PROTEIN NAME	Transcription factor p65
IMMUNOGEN	Synthesized peptide derived from human NFκB-p65 around the acetylation site of K310. at AA range: 250-330

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse
TESTED APPLICATIONS	WB,IHC,IF,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.