

Acetyl-eIF5A/eIF5A2 (K47) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10095	SYNONYMS
SIZE	100 ug	EIF5A2; Eukaryotic translation initiation factor 5A-2;
HOST	Rabbit	eIF-5A-2; eIF-5A2; Eukaryotic initiation factor 5A isoform 2

PRODUCT INFORMATION

GENE NAME	EIF5A/EIF5A2
PROTEIN NAME	Eukaryotic translation initiation factor 5A-2
IMMUNOGEN	Synthesized peptide derived from the human eIF5A/eIF5A2 around the acetylation site of K47.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.