

## Acetyl-HNF-3 $\alpha$ / $\beta$ / $\gamma$ (K264/253/211) Polyclonal Antibody

### OVERVIEW

---

<b>CATALOG NO.</b>	PTX10100	<b>SYNONYMS</b>
<b>SIZE</b>	100 ug	FOXA1; HNF3A; TCF3A; Hepatocyte nuclear factor 3-alpha; HNF-3-alpha; HNF-3A; Forkhead box protein A1; Transcription factor 3A; TCF-3A; FOXA2; HNF3B; TCF3B; Hepatocyte nuclear factor 3-beta; HNF-3-beta; HNF-3B; Forkhead box protein A2; Transcription factor
<b>HOST</b>	Rabbit	

### PRODUCT INFORMATION

---

<b>GENE NAME</b>	FOXA1/FOXA2/FOXA3
<b>PROTEIN NAME</b>	Hepatocyte nuclear factor 3-alpha/Hepatocyte nuclear factor 3-beta/Hepatocyte nuclear factor 3-gamma
<b>IMMUNOGEN</b>	Synthesized peptide derived from the human HNF-3 $\alpha$ / $\beta$ / $\gamma$ around the acetylation site of K264/253/211.

### APPLICATIONS AND STORAGE

---

<b>SPECIES REACTIVITY</b>	Human, Mouse, Rat
<b>TESTED APPLICATIONS</b>	WB,ELISA
<b>RECOMMENDED DILUTION</b>	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
<b>DELIVERY CONDITIONS</b>	Blue Ice
<b>STABILITY &amp; STORAGE</b>	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

**For research use only.**