

Acetyl-HSP 90 (K292/284) Polyclonal Antibody

OVERVIEW

| | | |
|--------------------|----------|---|
| CATALOG NO. | PTX10103 | SYNONYMS |
| SIZE | 100 ug | HSP90AA1; HSP90A; HSPC1; HSPCA; Heat shock protein HSP 90-alpha; Heat shock 86 kDa; HSP 86; HSP86; Renal carcinoma antigen NY-REN-38; |
| HOST | Rabbit | HSP90AB1; HSP90B; HSPC2; HSPCB; Heat shock protein HSP 90-beta; HSP 90; Heat shock 84 kDa; HSP 84; HSP84 |

PRODUCT INFORMATION

| | |
|---------------------|--|
| GENE NAME | HSP90AA1/HSP90AB1 |
| PROTEIN NAME | Heat shock protein HSP 90-alpha/Heat shock protein HSP 90-beta |
| IMMUNOGEN | Synthesized peptide derived from the human HSP 90 around the acetylation site of K292/284. |

APPLICATIONS AND STORAGE

| | |
|--------------------------------|--|
| SPECIES REACTIVITY | Human, Mouse, Rat |
| TESTED APPLICATIONS | WB,IHC,ELISA |
| RECOMMENDED DILUTION | Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications. |
| DELIVERY CONDITIONS | Blue Ice |
| STABILITY & STORAGE | PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles |

For research use only.