

Acetyl-HSP 90 (K435) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10104	SYNONYMS
SIZE	100 ug	HSP90AA1; HSP90A; HSPC1; HSPCA; Heat shock protein HSP 90-alpha; Heat shock 86 kDa; HSP 86;
HOST	Rabbit	HSP86; Renal carcinoma antigen NY-REN-38

PRODUCT INFORMATION

GENE NAME	HSP90AA1
PROTEIN NAME	Heat shock protein HSP 90-alpha
IMMUNOGEN	Synthesized peptide derived from the human HSP 90 around the acetylation site of K435.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.