

Acetyl-p73 (K327) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10112	SYNONYMS
SIZE	100 ug	TP73; P73; Tumor protein p73; p53-like transcription factor; p53-related protein
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	TP73
PROTEIN NAME	Tumor protein p73
IMMUNOGEN	Synthesized peptide derived from the human p73 around the acetylation site of K327.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.