

## Acetyl-Ub (K27) Polyclonal Antibody

### OVERVIEW

---

<b>CATALOG NO.</b>	PTX10122	<b>SYNONYMS</b>
<b>SIZE</b>	100 ug	UBB; Polyubiquitin-B; UBC; Polyubiquitin-C; RPS27A;
<b>HOST</b>	Rabbit	UBA80; UBCEP1; Ubiquitin-40S ribosomal protein S27a; Ubiquitin carboxyl extension protein 80; UBA52; UBCEP2; Ubiquitin-60S ribosomal protein L40; CEP52; Ubiquitin A-52 residue ribosomal protein fusion pr

### PRODUCT INFORMATION

---

<b>GENE NAME</b>	UBA52/RPS27A/UBB/UBC
<b>PROTEIN NAME</b>	Ubiquitin
<b>IMMUNOGEN</b>	Synthesized peptide derived from the human Ub around the acetylation site of K27.

### APPLICATIONS AND STORAGE

---

<b>SPECIES REACTIVITY</b>	Human, Mouse, Rat
<b>TESTED APPLICATIONS</b>	WB,IHC,ELISA
<b>RECOMMENDED DILUTION</b>	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
<b>DELIVERY CONDITIONS</b>	Blue Ice
<b>STABILITY &amp; STORAGE</b>	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

**For research use only.**