

Acetyl-PCAF (K428) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10126	SYNONYMS
SIZE	100 ug	KAT2B; PCAF; Histone acetyltransferase KAT2B;
HOST	Rabbit	Histone acetyltransferase PCAF; Histone acetylase PCAF; Lysine acetyltransferase 2B; P300/CBP- associated factor; P/CAF

PRODUCT INFORMATION

GENE NAME	KAT2B
PROTEIN NAME	Histone acetyltransferase KAT2B
IMMUNOGEN	Synthesized peptide derived from the human PCAF around the acetylation site of K428.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.