

Acetyl-PTEN (K402) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10127	SYNONYMS
SIZE	100 ug	PTEN; MMAC1; TEP1; Phosphatidylinositol 3,4,5-
HOST	Rabbit	trisphosphate 3-phosphatase and dual-specificity protein phosphatase PTEN; Mutated in multiple advanced cancers 1; Phosphatase and tensin homolog

PRODUCT INFORMATION

GENE NAME	PTEN
PROTEIN NAME	Phosphatidylinositol 3,4,5-trisphosphate 3-phosphatase and dual-specificity protein phosphatase PTEN
IMMUNOGEN	Synthesized peptide derived from the human PTEN around the acetylation site of K402.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.