

Acetyl-Ub (K29) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10128	SYNONYMS
SIZE	100 ug	UBB; Polyubiquitin-B; UBC; Polyubiquitin-C; RPS27A;
HOST	Rabbit	UBA80; UBCEP1; Ubiquitin-40S ribosomal protein S27a; Ubiquitin carboxyl extension protein 80; UBA52; UBCEP2; Ubiquitin-60S ribosomal protein L40; CEP52; Ubiquitin A-52 residue ribosomal protein fusion pr

PRODUCT INFORMATION

GENE NAME	UBA52/RPS27A/UBB/UBC
PROTEIN NAME	Ubiquitin
IMMUNOGEN	Synthesized peptide derived from the human Ub around the acetylation site of K29.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.