

Acetyl-CRM1 (K568) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10132	SYNONYMS
SIZE	100 ug	XPO1; CRM1; Exportin-1; Exp1; Chromosome region maintenance 1 protein homolog
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	XPO1
PROTEIN NAME	Exportin-1
IMMUNOGEN	Synthesized peptide derived from the Internal region of human CRM1 around the acetylation site of K568.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.