

β -actin Rabbit Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX1024
SIZE	100 μ L
HOST	Rabbit
ISOTYPE	IgG
CLONALITY	Polyclonal Antibody

APPLICATIONS AND STORAGE

APPLICATIONS	WB, IHC, IF
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles
DELIVERY CONDITION	Blue ice