

β II tubulin Monoclonal Antibody

OVERVIEW

CATALOG NO.	PTX10290	SYNONYMS
SIZE	100 ug	
HOST	Mouse	

PRODUCT INFORMATION

GENE NAME	TUBB2A
PROTEIN NAME	Tubulin beta-2A chain
IMMUNOGEN	Synthetic Peptide

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Rat, (Human)
TESTED APPLICATIONS	WB, IP
RECOMMENDED DILUTION	WB: 1:50000-100000 IP:1:200
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.