

FH Fumarase Monoclonal Antibody

OVERVIEW

CATALOG NO.	PTX10295	SYNONYMS
SIZE	100 ug	
HOST	Mouse	

PRODUCT INFORMATION

GENE NAME	FH
PROTEIN NAME	Fumarate hydratase, mitochondrial
IMMUNOGEN	Recombinant Protein

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB, IF
RECOMMENDED DILUTION	WB: 1:1000-2000
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.