

Flotillin-1 Monoclonal Antibody

OVERVIEW

CATALOG NO.	PTX10296	SYNONYMS
SIZE	100 ug	FN1; FN; Fibronectin; FN; Cold-insoluble globulin; CIG
HOST	Mouse	

PRODUCT INFORMATION

GENE NAME	FN1
PROTEIN NAME	Fibronectin
IMMUNOGEN	Synthetic Peptide

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Rat, (Human)
TESTED APPLICATIONS	WB
RECOMMENDED DILUTION	WB: 1:1000-2000
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.