

COX2 Monoclonal Antibody

OVERVIEW

CATALOG NO.	PTX10298	SYNONYMS
SIZE	100 ug	PTGS2; COX2; Prostaglandin G/H synthase 2;
HOST	Mouse	Cyclooxygenase-2; COX-2; PHS II; Prostaglandin H2 synthase 2; PGH synthase 2; PGHS-2; Prostaglandin-endoperoxide synthase 2

PRODUCT INFORMATION

GENE NAME	PTGS2
PROTEIN NAME	Prostaglandin G/H synthase 2
IMMUNOGEN	Synthetic Peptide

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	IHC
RECOMMENDED DILUTION	IHC: 1:200
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.