

## Histone H2B Monoclonal Antibody

### OVERVIEW

---

<b>CATALOG NO.</b>	PTX10301	<b>SYNONYMS</b>
<b>SIZE</b>	100 ug	HIST1H2BB; H2BFF; Histone H2B type 1-B; Histone
<b>HOST</b>	Mouse	H2B.1; Histone H2B.f; H2B/f

### PRODUCT INFORMATION

---

<b>GENE NAME</b>	HIST1H2BB
<b>PROTEIN NAME</b>	Histone H2B type 1-B
<b>IMMUNOGEN</b>	Synthetic Peptide

### APPLICATIONS AND STORAGE

---

<b>SPECIES REACTIVITY</b>	Human, Mouse, Rat
<b>TESTED APPLICATIONS</b>	WB
<b>RECOMMENDED DILUTION</b>	WB: 1:1000-3000
<b>DELIVERY CONDITIONS</b>	Blue Ice
<b>STABILITY &amp; STORAGE</b>	PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

**For research use only.**