

eNOS Monoclonal Antibody

OVERVIEW

CATALOG NO.	PTX10303	SYNONYMS
SIZE	100 ug	NOS3; Nitric oxide synthase, endothelial; Constitutive
HOST	Mouse	NOS; cNOS; EC-NOS; Endothelial NOS; eNOS; NOS type III; NOSIII

PRODUCT INFORMATION

GENE NAME	NOS3
PROTEIN NAME	Nitric oxide synthase, endothelial
IMMUNOGEN	Recombinant Protein

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB
RECOMMENDED DILUTION	WB: 1:500-2000
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.