

## GAPDH Polyclonal Antibody

### OVERVIEW

---

<b>CATALOG NO.</b>	PTX10339	<b>SYNONYMS</b>
<b>SIZE</b>	100 ug	GAPDH; GAPD; CDABP0047; OK/SW-cl.12;
<b>HOST</b>	Rabbit	Glyceraldehyde-3-phosphate dehydrogenase; GAPDH; Peptidyl-cysteine S-nitrosylase GAPDH

### PRODUCT INFORMATION

---

<b>GENE NAME</b>	GAPDH
<b>PROTEIN NAME</b>	Glyceraldehyde-3-phosphate dehydrogenase
<b>IMMUNOGEN</b>	Recombinant Protein

### APPLICATIONS AND STORAGE

---

<b>SPECIES REACTIVITY</b>	Zebrafish
<b>TESTED APPLICATIONS</b>	WB
<b>RECOMMENDED DILUTION</b>	WB: 1:2000-5000
<b>DELIVERY CONDITIONS</b>	Blue Ice
<b>STABILITY &amp; STORAGE</b>	PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

**For research use only.**