

## **Bovine Serum Albumin Polyclonal Antibody**

### **OVERVIEW**

---

<b>CATALOG NO.</b>	PTX10365	<b>SYNONYMS</b>
<b>SIZE</b>	100 ug	
<b>HOST</b>	Rabbit	

### **PRODUCT INFORMATION**

---

**GENE NAME**

**PROTEIN NAME**

**IMMUNOGEN**      Recombinant Protein

### **APPLICATIONS AND STORAGE**

---

<b>SPECIES REACTIVITY</b>	Cow
<b>TESTED APPLICATIONS</b>	WB, ELISA
<b>RECOMMENDED DILUTION</b>	WB: 1:2000-5000
<b>DELIVERY CONDITIONS</b>	Blue Ice
<b>STABILITY &amp; STORAGE</b>	PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

**For research use only.**