

Histone H3 (Acetyl Lys4) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10379	SYNONYMS
SIZE	100 ug	HIST1H3A; H3FA; HIST1H3B; H3FL; HIST1H3C; H3FC; HIST1H3D; H3FB; HIST1H3E; H3FD; HIST1H3F; H3FI; HIST1H3G; H3FH; HIST1H3H; H3FK; HIST1H3I; H3FF; HIST1H3J; H3FJ; Histone H3.1; Histone H3/a; Histone H3/b; Histone H3/c; Histone H3/d; Histone H3/f; Histone H3/
HOST	Rabbit	

PRODUCT INFORMATION

GENE	HIST1H3A/HIST1H3B/HIST1H3C/HIST1H3D/HIST1H3E/HIST1H3F/HIST1H3G/HIST1H3H/HIST1H3I/ HIST1H3J/HIST2H3A/HIST2H3C/HIST2H3D/H3F3A/H3F3B
PROTEI	Histone H3.1/Histone H3.2/Histone H3.3
N NAME	
IMMUN	Synthetic Peptide
OGEN	

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB
RECOMMENDED DILUTION	WB: 1:500-1000
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.