

Histone H2B (Acetyl Lys15) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10413	SYNONYMS
SIZE	100 ug	HIST1H2BA; TSH2B; Histone H2B type 1-A; Histone H2B, testis; Testis-specific histone H2B; HIST1H2BB; H2BFF; Histone H2B type 1-B; Histone H2B.1; Histone H2B.f; H2B/f; HIST1H2BC; H2BFL; HIST1H2BE; H2BFH; HIST1H2BF; H2BFG; HIST1H2BG; H2BFA; HIST1H2BI; H2BFK
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	HIST1H2BC
PROTEIN NAME	Histone H2B type 1-A/Histone H2B type 1-B/Histone H2B type 1-C/E/F/G/I
IMMUNOGEN	Synthetic Peptide

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB
RECOMMENDED DILUTION	WB: 1:200-500
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.