

Phospho-Cdc2 (Y15) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10463	SYNONYMS
SIZE	100 ug	CDK1; CDC2; CDC28A; CDKN1; P34CDC2; Cyclin-dependent kinase 1; CDK1; Cell division control protein 2 homolog; Cell division protein kinase 1; p34 protein kinase
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	CDK1
PROTEIN NAME	Cyclin-dependent kinase 1
IMMUNOGEN	Synthesized peptide derived from human Cdc2 around the phosphorylation site of Y15.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat, Monkey
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.