

Phospho-Chk2 (S516) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10469	SYNONYMS
SIZE	100 ug	CHEK2; CDS1; CHK2; RAD53; Serine/threonine-
HOST	Rabbit	protein kinase Chk2; CHK2 checkpoint homolog; Cds1 homolog; Hucds1; hCds1; Checkpoint kinase 2

PRODUCT INFORMATION

GENE NAME	CHEK2
PROTEIN NAME	Serine/threonine-protein kinase Chk2
IMMUNOGEN	Synthesized peptide derived from human Chk2 around the phosphorylation site of S516.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Monkey
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.