

Phospho-DARPP-32 (T75) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10479	SYNONYMS
SIZE	100 ug	PPP1R1B; DARPP32; Protein phosphatase 1 regulatory subunit 1B; DARPP-32; Dopamine- and cAMP-regulated neuronal phosphoprotein
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	PPP1R1B
PROTEIN NAME	Protein phosphatase 1 regulatory subunit 1B
IMMUNOGEN	Synthesized peptide derived from human DARPP-32 around the phosphorylation site of T75.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat, Monkey
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.