

Phospho-EGFR (Y1092) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10483	SYNONYMS
SIZE	100 ug	EGFR; ERBB; ERBB1; HER1; Epidermal growth factor
HOST	Rabbit	receptor; Proto-oncogene c-ErbB-1; Receptor tyrosine-protein kinase erbB-1

PRODUCT INFORMATION

GENE NAME	EGFR
PROTEIN NAME	Epidermal growth factor receptor
IMMUNOGEN	Synthesized peptide derived from human EGFR around the phosphorylation site of Y1092.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.