

Phospho-FAK (Y861) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10494	SYNONYMS
SIZE	100 ug	PTK2; FAK; FAK1; Focal adhesion kinase 1; FADK 1;
HOST	Rabbit	Focal adhesion kinase-related nonkinase; FRNK; Protein phosphatase 1 regulatory subunit 71; PPP1R71; Protein-tyrosine kinase 2; p125FAK; pp125FAK

PRODUCT INFORMATION

GENE NAME	PTK2
PROTEIN NAME	Focal adhesion kinase 1
IMMUNOGEN	Synthesized peptide derived from human FAK around the phosphorylation site of Y861.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat, Monkey
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.