

Phospho-FoxO4 (S197) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10498	SYNONYMS
SIZE	100 ug	FOXO4; AFX; AFX1; MLLT7; Forkhead box protein
HOST	Rabbit	O4; Fork head domain transcription factor AFX1

PRODUCT INFORMATION

GENE NAME	FOXO4
PROTEIN NAME	Forkhead box protein O4
IMMUNOGEN	Synthesized peptide derived from human FoxO4 around the phosphorylation site of S197.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.