

Phospho-GATA-1 (S310) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10501	SYNONYMS
SIZE	100 ug	GATA1; ERYF1; GF1; Erythroid transcription factor;
HOST	Rabbit	Eryf1; GATA-binding factor 1; GATA-1; GF-1; NF-E1 DNA-binding protein

PRODUCT INFORMATION

GENE NAME	GATA1
PROTEIN NAME	Erythroid transcription factor
IMMUNOGEN	Synthesized peptide derived from human GATA-1 around the phosphorylation site of S310.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat, Monkey
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.