

Phospho-HDAC8 (S39) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX10507	SYNONYMS
SIZE	100 ug	HDAC8; HDAC1; CDA07; Histone deacetylase 8;
HOST	Rabbit	HD8

PRODUCT INFORMATION

GENE NAME	HDAC8
PROTEIN NAME	Histone deacetylase 8
IMMUNOGEN	Synthesized peptide derived from human HDAC8 around the phosphorylation site of S39.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.